

PRODUCT SPECIFICATION

Product Name : NuHospisol-G – Disinfecting and sterilizing solution with Activator

Special Features

- For effective, safe rapid disinfection and sterilization of endoscopes, surgical instruments, anaesthetic equipment, urological instruments, neonatal care instruments, dental instruments etc.
- Contains added rust inhibitors

Applications

- Instrument disinfectant for cold sterilization where autoclaving is not feasible
- Especially recommended for use with : lensed instruments (flexible fiber-optic and rigid endoscopes) Metal, glad, rubber and plastic objects dental instruments

Product Information

Composition : Glutaraldehyde solution 2% w/v
Labeling : Shall comply with the corresponding labeling specification and customer requirement.
Shelf Life : 3 years for concentrated & 14 days for activated solution.
Storage condition : Shall be stored in cool dry place and away from direct light, away from acids, oxidizing agents & heat.

Packaging Details

Packing	Units per Carton
1 Litre Can with Activator Vial	12 Nos
5 Litre Can with Activator Vial	2 Nos.

Certifications : ISO 9001 : 2008
ISO 13485 : 2003

Usage & Dilution

To be used undiluted after activation Activated solution may be reused

Activity Time

Standard Disinfection	10 Minutes
High Level Disinfection	45 Minutes
Sterilization	10 Hours

Antimicrobial Activity

Bacteria	Yes
Fungus	Yes
Virus	Yes
TB	Yes
Spores	Yes