

PRODUCT SPECIFICATION

PRODUCT NAME : VINYL EXAMINATION GLOVES NON STERILE POWDERFREE

Special Features

- Latex-free, odour-free and powder-free
- Can withstand IPA washing (ISO Propyl Alcohol).
- Can be used in food industry
- Ideal for doctors with latex allergy

Product Information

Material : 100% Synthetic Poly Vinyl Chloride (PVC), No Natural Rubber Latex
 Colour : Clear
 Design : Ambidextrous, non-sterile, beaded cuff and smooth external finish.
 Powder Content : ≤ 2.0 mg max per one piece of gloves
 Sterilization : Non Sterile
 Labeling : Shall comply with the corresponding labeling specification and customer requirement.
 Shelf Life : 3 years from the date of manufacture
 Storage condition : Shall be stored in cool dry place and away from direct light

Packaging Details

Inner Box	100 pcs
Master Carton	1000 pcs (10 inner boxes)
20 feet container	15,00,000 pcs (1500 cartons)
40 feet HC container	35,00,000 pcs (3500 cartons)

Certifications : ISO 9001 : 2008
 ISO 13485 : 2003

Physical Dimensions

Length(mm) Min	Size	XS	S	M	L	XL
	Length	Min 230mm				
Width (mm)	Size	XS	S	M	L	XL
	Width (±6)	75mm	85mm	95mm	105mm	115mm
Weight	(+/- 0.5 gm)	3.7gm	3.8gm	4.2gm	4.4gm	4.7gm
Thickness(mm) Min	Cuff	--				
	Palm	Min 0.05mm (Single Wall)				
	Finger	Min 0.05mm (Single Wall)				

Physical Properties

Characteristics	Before Ageing	After Ageing (70± 2 °C for 166 ± 2 hrs)
Tensile Strength (Mpa) min.	9 mpa Min	9 mpa Min
Ultimate Elongation (%) min.	300% Min	300% Min

Inspection Level (Pre-Shipment)

Freedom from Holes : Sampling as per ISO 2859, Inspection level - I & AQL - 1.5 or AQL - 2.5
(ASTM D 5151)

Visual defects-Major : Inspection level - I & AQL : 1.5

Visual defects-Minor : Inspection level – I & AQL: 2.5

Directions for use:

1.1 Preparation For Use : Ensure selection of gloves with correct size.

Below is the correct glove size width:

Size	XS	S	M	L	XL
Width Range	75mm	85mm	95mm	105mm	115mm

1.2 How To Use

- Dry hands before donning the gloves.
- Take out the product after checking package.
- Gloves need to be donned carefully; do not apply excessive force to pull the gloves when donning.
- Please use it in time after taking out from the box, and avoid contacting with ketones.
- Keep the fingernails short and change gloves frequently to avoid tearing issue.
- The gloves should be discarded with an appropriate waste disposal container after use according to local disposal regulations.

1.3 Cautions For Use

- The product is for single use only.
- Does not contain natural rubber latex.

1.4 Storage Condition

- The atmosphere must always be conducive to the storage of the gloves.
- Gloves before and after opened for use shall be stored in dark, cool and dry place.
- Gloves need to avoid area with excessive heat. The storage areas need to be controlled below 40°C/ 104 °F.
- Gloves need to be stored at the areas with good air ventilation.
- Gloves need to be shielded from direct sunlight and fluorescent lighting.
- Gloves need to be stored far away from any energy radiating and heat generating equipment.

1.5 Expiration Date

The shelf life of the gloves is 3 years from the date of manufacture.