

Nulife

Nulife is a global medical disposable solutions company with its footprint in 80+ countries and its headquarters in Mumbai, India.

With a workforce of over 2000 people worldwide, Nulife annually sells more than 300 million medical devices.

Triple Safety

Nulife Triple Safety is a brand that is committed to offering affordable healthcare solutions.

Powered by cutting-edge innovation, Nulife Triple Safety has disrupted the status quo in its sector through cost-effective healthcare products over the years.



Nulife Triple Safety Latex Examination Gloves (Pre-powdered or powder-free)

Product Information

Nulife Triple Safety Latex Examination Gloves are designed for multipurpose use.

These gloves are available in pre-powdered as well as powder-free options.



They have been tested and proven safe against dermal sensitisation as well as skin irritation.



Special Features

- Diverse applications in medical and therapeutic procedures
- Suitable for handling infectious medical materials
- Available in pre-powdered or powder-free options

Latex Gloves: Specifications

Latex Examination Gloves, Pre-powdered or Powder-free, Non-sterile						Standard Reference
Size	Extra Small	Small	Medium	Large	Extra Large	EN 455-2: 2015
Length (mm)	≥ 240					EN 455-2: 2015
Width (mm)	≤ 80	80 ± 10	95 ± 10	110 ± 10	≥ 110	EN 455-2: 2015
Tensile Force (N)	≥ 06					EN 455-2: 2015
Powder Content	Powder-free (Max. 2mg/glove) and Powdered (Max. 10mg/dm ²)					EN 455-3: 2015 & ASTM 3578-19 & ASTM D 6124-06
Protein Leachable	Max. 200µg/dm ²					EN 455-3:2015 ASTM 3578-19
Freedom from holes	Should pass leak test					EN 455-1: 2020

Disclaimer:

1. The description and specifications on the website are correct at the time of publishing.
2. In pursuit of continual improvement, we reserve the right to change the design and specifications of any of our products without prior notice.
3. All photographs, product designs and trademarks represented in this website are the sole intellectual property of MRK HEALTHCARE PVT. LTD.