

Nulife

Nulife is a global medical disposable solutions company with its footprint in 80+ countries and its headquarters in Mumbai, India.

With a workforce of over 2000 people worldwide, Nulife annually sells more than 300 million medical devices.

Triple Safety

Nulife Triple Safety is a brand that is committed to offering affordable healthcare solutions.

Powered by cutting-edge innovation, Nulife Triple Safety has disrupted the status quo in its sector through cost-effective healthcare products over the years.



Nulife Triple Safety Nitrile Examination Gloves

(Powder-free, Latex-free)

Product Information

Nulife Triple Safety Nitrile Examination Gloves are latex-free and are manufactured using synthetic rubber and are ideal for people with latex allergy.



The gloves are solvent-resistant and odour-free, making them useful in the food and dairy industry.



Special Features

- Prevents allergic user reactions
- Ambidextrous usage
- Provides better resistance to chemicals than NRL gloves

Nitrile Gloves: Specifications

Nitrile Examination Gloves (Latex-free), Powder-free, Non-sterile						Standard Reference
Size	Extra Small	Small	Medium	Large	Extra Large	
Length (mm)	≥ 240					EN 455-2: 2015
Width (mm)	≤ 80	80 ± 10	95 ± 10	110 ± 10	≥ 110	EN 455-2: 2015
Tensile Force (N) before & after ageing	≥ 06					EN 455-2: 2015
Powder Content	Max. 2mg/glove					EN 455-3: 2015 & ASTM 3578-19 & ASTM D 6124-06
Freedom from holes	Should pass leak test					EN 455-1: 2020

Disclaimer:

1. The description and specifications on the website are correct at the time of publishing.
2. In pursuit of continual improvement, we reserve the right to change the design and specifications of any of our products without prior notice.
3. All photographs, product designs and trademarks represented in this website are the sole intellectual property of MRK HEALTHCARE PVT. LTD.