

Nulife

Nulife is a global medical disposable solutions company with its footprint in 80+ countries and its headquarters in Mumbai, India.

With a workforce of over 2000 people worldwide, Nulife annually sells more than 300 million medical devices.

Triple Safety

Nulife Triple Safety is a brand that is committed to offering affordable healthcare solutions.

Powered by cutting-edge innovation, Nulife Triple Safety has disrupted the status quo in its sector through cost-effective healthcare products over the years.



Nulife Triple Safety Sterile Latex Surgical Gloves (Powder-free)

Product Information

Nulife Triple Safety Sterile Latex Surgical Gloves (Powder-free) offer exceptional dexterity, enhance comfort and flexibility.

The gloves have minimal extractable NR latex protein content, minimising the risk of skin allergies.

Furthermore, their micro-rough texture ensures a secure grip.



Special Features

- Soft and flexible
- Unique HiLeach process minimises protein content leading to lower allergies
- Optimum moisture content

TECHNICAL SPECIFICATION SHEET

Product Information

Product Name	Nulife Triple Safety Sterile Latex Surgical Gloves (Powder-free)
Material	Natural Rubber Latex
Colour	Creamy White
Design	Hand-specific, micro-rough textured
Type	Powder-free Latex Surgical Gloves
Sterilisation	Ethylene Oxide, Gamma upon request
Labelling	Shall comply with the corresponding labelling specification, customer requirement and ISO EN15223-1: 2021
Shelf Life	5 years from the date of manufacture
Storage Condition	Store in cool, clean & dry condition, away from heat & direct sunlight. Recommended storage condition: 5 °C to 40 °C
Certifications	ISO 13485: 2016, CE (MDR (EU 2017/745), PPE (EU 2016/425), USFDA(510 K)
Shape	Anatomic
Cuff	Beaded
External Surface	Palm-textured
Internal Surface	Powder-free

Physical Dimension

Size	Minimum Length (mm)		Palm Width (mm)		Minimum Thickness (mm)		
	EN 455-2	MRK	EN 455-2	MRK	Location	ASTM D 3577	MRK
5.5	250	280 Min	72±4	71-76	Cuff	0.1	0.11 Min
6	260	280 Min	77±5	77-82			
6.5	260	280 Min	83±5	82-86	Palm	0.1	0.15 Min
7	270	280 Min	89±5	88-94			
7.5	270	280 Min	95±5	95-99	Finger	0.1	0.16 Min
8	270	280 Min	102±6	101-107			
8.5	280	280 Min	108±6	108-111			
9	280	280 Min	114±6	114-115			

Physical Properties

Characteristics	Before Ageing		After Ageing (100±2 °C for 22±0.3 h/70 °C for 7 days)	
	ASTM D 3577	MRK	ASTM D 3577	MRK
Tensile Strength (MPa)	24 min	24 min	18 min	18 min
Ultimate Elongation (%)	750 min	760 min	560 min	600 min
Stress at 500% Elongation (MPa)	5.5 max	5.0 max	N/A	N/A
Force at Break (N)	EN 455-2	MRK	EN 455-2	MRK
	9.0 min	9.0 min	9.0 min	9.0 min

Freedom from holes	Inspection Level G1 AQL 0.65
Protein Content	Max. 200µg/dm ² as per ASTM D 3577
Powder Content	Max. 2 mg/glove
Packaging	10x50 = 500 pairs
Shelf Life Period	5 years
CE-classification	MDR: Class IIa, PPE: Category III
Compliances	EN 455 parts 1,2,3 & 4, ASTM D 3577, ISO 10282, IS 13422, IS 4148, ISO 13485: 2016, CE (MDR (EU 2017/745), USFDA(510 K)
Compliances PPE	PPE (EU 2016/425), EN ISO 374-1,2,4,5, EN 16523-1, EN ISO 21420
Notified Body	MDR: DNV Product Assurance AS (2460), PPE: SGS Fimko oy (NB 0598)

Disclaimer:

1. The description and specifications on the website are correct at the time of publishing.
2. In pursuit of continual improvement, we reserve the right to change the design and specifications of any of our products without prior notice.
3. All photographs, product designs and trademarks represented in this website are the sole intellectual property of MRK HEALTHCARE PVT. LTD.

MRK HEALTHCARE PVT LTD

A 1201, Naman Midtown, Senapati Bapat Marg,
 Elphinstone Road, Mumbai 400 013, Maharashtra. INDIA
 T: +91 22 6261 1111 | F: +91 22 6261 1199
 E: info@mrkhealthcare.com